



8006

Aneroid sphygmomanometer



FEATURES

Sphygmomanometers with nylon cuff with velcro closure, reliable manometer with clip to be fixed on the cuff. Scale 0-300 mmHg. Floating zero, nylon velcro cuff, latex bulb with chromium plated valve. Chromium-plated air tap with air release regulation by screw valve. Supplied with a simulated-leather soft pouch with zip.

MAIN FEATURES

Size

Adults

SEE ALSO



8004
Colored aneroid
sphygmomanometer with carrying
case



8012
Red palm-style aneroid
sphygmomanometer plus



1740
Rotating wall aneroid
sphygmomanometer

ACCESSORIES



8143
Double tube pediatric cuff



8141
Adults cuff with double tube



8142
Single tube pediatric cuff



8140
Adults cuff with single tube

ACCESSORY OF



1366
Red knapsack "115" with side
pockets and compartments
COMPLETE
