



8008

One hand aneroid sphygmomanometer



FEATURES

Antishock ABS case manometer \varnothing 60 mm. Range 0-300 mm/Hg. Accuracy \pm 3 mm/Hg. Manometer with NON STOP PIN . Easy left-hand or right-hand operation. Adult velcro cuff with metal ring. Standard model in vinyl case.

MAIN FEATURES

Size

Adults

SEE ALSO



8012
Red palm-style aneroid
sphygmomanometer plus



1740
Rotating wall aneroid
sphygmomanometer



8004
Colored aneroid
sphygmomanometer with carrying
case



8006
Aneroid sphygmomanometer

ACCESSORIES



8142
Single tube pediatric cuff



8141
Adults cuff with double tube



8143
Double tube pediatric cuff



8140
Adults cuff with single tube
