



8008

One hand aneroid sphygmomanometer



FEATURES

Classic, lightweight, reliable palm model aneroid sphygmo with certified manometer ABS housing, integral PVC bulb and air release valve for one-hand use, provided with handle, blue nylon, velcro closure, calibrated cuff and 1 tube TPU bladder.

Easy to read: large dial \emptyset 60 mm, with rubber protective ring.

Supplied in a nice nylon carrying bag with external velcro pocket and cardboard box. Latex free.

MAIN FEATURES			
Size	Adults		

SEE ALSO



8012Red palm-style aneroid sphygmomanometer plus



8004
Colored aneroid sphygmomanometer with carrying case



ACCESSORIES







