



8014

Compact aneroid sphigmomanometer



FEATURES

The handmade sphigmomanometer has high quality finish. Its components are made applying the highest technological solutions: German movement follows the highest quality standards. Polycarbonate dial glass, alcohol and mechganical strain proof. Sphigmomanometer drain valve easy to activate and very sensitive. Natural rubber bladder. Cuff made by a strong anti-allergic cotton. Manometer, fixed on the cuff, is highly resistant to falls or casual blows. Made in Italy according to the EN 1060-1 and EN 1060-2 rules.

M.	ΔΙ	N	FF	ΞΔ٦	ΓIJ	DI	ES
141	\sim 1	I N	гь		ı	T.	

Size Adults

SEE ALSO



8008 Solved aneroid Sphygmomanometer Sph