



## 8020/G

---

Yellow flat phonendoscope



---

### FEATURES

Flat phonendoscope with anodized stamped and lathed aluminium head diameter 46 mm. Threaded metal ring for membrane blocking. Metal chromium plated binaural, length 150 mm. Y rubber tube in one piece, diameter 5 mm length 60 cm. Single printed carton box. Available in 3 different colors.

---

### MAIN FEATURES

Size	Adults
Cardio	

---

SEE ALSO

---



**8030/B**  
Blue stethophonendoscope  
"Rappaport"



**8025/N**  
Black stethophonendoscope



**8020/N**  
Black flat phonendoscope



**8025/G**  
Yellow stethophonendoscope



**8020/R**  
Red flat phonendoscope



**8025/R**  
Red stethophonendoscope

---