



# 8025/N

Black stethophonendoscope



## FEATURES

Stethophonendoscope with head diameter 47/31 mm. Aluminium anodized head, chromium plated metal binaural, Y rubber tube in one piece. Rotating piper-fitting for both functions. Available in 3 different colors.

## MAIN FEATURES

Size	Adults
Cardio	<input checked="" type="checkbox"/>

## SEE ALSO



**8020/G**  
Yellow flat phonendoscope



**8025/G**  
Yellow stethophonendoscope



**8030/B**  
Blue stethophonendoscope  
"Rappaport"



**8025/R**  
Red stethophonendoscope



**8020/N**  
Black flat phonendoscope



**8020/R**  
Red flat phonendoscope