



## 8284

---

Multiparameter monitor MP1000



---

### FEATURES

---

Maximum compactness and portability thanks to reduced dimensions and weight

- Handle for a safe grip and an easy transport
- Easy usage by mean of pop-up menu
- ECG – SpO2 – NIBP – TEMP. monitoring
- Analysis of arrhythmias ST and pacemaker detection
- Allarm set up for all parameters
- 72 hours trend of cardiac freq., SpO2 e NISBP
- Operation from main or integrated rechargeable batteries
- Integrated thermic printer

**Standard accessories:** ECG cable (7 derivations)• Re-usable adult SpO2 sensor• Re-usable adult / pediatric NIBP cuff• NIBP tube• Operator manual **Optional accessories:** Temperature probe