



# 836 Patriot

Adjustable cervical collar



## FEATURES

Cervical collar adjustable to 4 different sizes. Designed to have in one article different possibilities of immobilization, it has been developed for the emergency medical use. Pre-molded chin support, locking clips, rear ventilation panel, enlarged trachea opening. The Patriot cervical collar is produced with high density polyethylene and foam padding; the flat one-piece design enables efficient storage where space is limited. X-ray lucent and easy to clean and disinfect.

## MAIN FEATURES

Size	Adults
Usage	Immobilization
Trauma	<input checked="" type="checkbox"/>

## SEE ALSO



**Patriot Baby 837**  
Adjustable cervical collar pediatric