



840 Antares

Rigid cervical collar in 2 pieces size INFANT



FEATURES

ANTARES cervical collar is a highly professional instrument indispensable for the emergency operator. It is composed of 2 pieces to be put on the front and back part of the neck enabling thus the perfect immobilization of the rachis. Central hole for tracheotomy and easy positioning thanks to the rip-off fastening which different colours indicate the various sizes. It is made with transparent polyethylene with E.V.A. washable inner lining fixed by plastic rivets. Available in 6 different sizes.

MAIN FEATURES

Size	Children
Usage	Immobilization
Trauma	
Pediatric	
Length (cm)	24,5
Max length (cm)	25
Max height (cm)	17
Min height (cm)	14,5

SEE ALSO



Antares 854

Cervical collar set 4 pieces with carrying case Art. 842 - 844 - 846 - 848



Sirius 876

Cervical collar set 6 pieces with carrying case