



## 842 Antares

Rigid cervical collar in 2 pieces size PEDIATRIC



### FEATURES

ANTARES cervical collar is an highly professional instrument indispensable for the emergency operator. It is composed by 2 pieces to be put on the front and back part of the neck enabling thus the perfect immobilization of the rachis. Central hole for tracheotomy and easy positioning thanks to the rip-off fastening which different colours indicate the various sizes. It is made with transparent polyethylene with E.V.A. washable inner lining fixed by plastic rivets. Available in 6 different sizes.

### MAIN FEATURES

Size	Children
Usage	Immobilization
Trauma	
Pediatric	
Max length (cm)	25
Length (cm)	24,5
Min height (cm)	13,5
Max height (cm)	14

### SEE ALSO



**Antares 856**  
Cervical collar set 6 pieces with carrying case



**Antares 854**  
Cervical collar set 4 pieces with carrying case Art. 842 - 844 - 846 - 848



**Sirius 874**  
Cervical collar set 4 pieces with carrying case Art. 862 - 864 - 866 - 868



**Sirius 876**  
Cervical collar set 6 pieces with carrying case