



## 846 Antares


Rigid cervical collar in 2 pieces size STANDARD



### FEATURES

ANTARES cervical collar is an highly professional instrument indispensable for the emergency operator. It is composed by 2 pieces to be put on the front and back part of the neck enabling thus the perfect immobilization of the rachis. Central hole for tracheotomy and easy positioning thanks to the rip-off fastening which different colours indicate the various sizes. It is made with transparent polyethylene with E.V.A. washable inner lining fixed by plastic rivets. Available in 6 different sizes.

### MAIN FEATURES

Size	Adults
Usage	Immobilization
Trauma	
Length (cm)	29
Max height (cm)	19
Min height (cm)	15

### SEE ALSO



#### Sirius 874

Cervical collar set 4 pieces with carrying case Art. 862 - 864 - 866 - 868



#### Sirius 876

Cervical collar set 6 pieces with carrying case



#### Antares 856

Cervical collar set 6 pieces with carrying case



#### Antares 854

Cervical collar set 4 pieces with carrying case Art. 842 - 844 - 846 - 848



#### Antares 854

Cervical collar set 4 pieces with carrying case Art. 842 - 844 - 846 - 848



#### Antares 854

Cervical collar set 4 pieces with carrying case Art. 842 - 844 - 846 - 848