



850 Antares

Rigid cervical collar in 2 pieces size TALL



FEATURES

ANTARES cervical collar is a highly professional instrument indispensable for the emergency operator. It is composed of 2 pieces to be put on the front and back part of the neck enabling thus the perfect immobilization of the rachis. Central hole for tracheotomy and easy positioning thanks to the rip-off fastening which different colours indicate the various sizes. It is made with transparent polyethylene with E.V.A. washable inner lining fixed by plastic rivets. Available in 6 different sizes.

MAIN FEATURES

Size	Adults
Usage	Immobilization
Trauma	
Length (cm)	35
Min height (cm)	20
Max height (cm)	22,5

SEE ALSO



Antares 854

Cervical collar set 4 pieces with carrying case Art. 842 - 844 - 846 - 848



Sirius 876

Cervical collar set 6 pieces with carrying case



Sirius 874

Cervical collar set 4 pieces with carrying case Art. 862 - 864 - 866 - 868



Antares 856

Cervical collar set 6 pieces with carrying case