



## 854 Antares

---

Cervical collar set 4 pieces with carrying case

Art. 842 - 844 - 846 - 848



---

### FEATURES

ANTARES cervical collar is a highly professional instrument indispensable for the emergency operator. It is composed of 2 pieces to be put on the front and back part of the neck enabling thus the perfect immobilization of the rachis. Central hole for tracheotomy and easy positioning thanks to the rip-off fastening which different colours indicate the various sizes. It is made with transparent polyethylene with E.V.A. washable inner lining fixed by plastic rivets. Available in 6 different sizes.

Art. 842 - 844 - 846 - 848

---

### MAIN FEATURES

Usage

Trauma

Immobilization



---

SEE ALSO

---



**Sirius 874**  
Cervical collar set 4 pieces with  
carrying case Art. 862 - 864 - 866 -  
868



**Sirius 876**  
Cervical collar set 6 pieces with  
carrying case



**Antares 844**  
Rigid cervical collar in 2 pieces size  
STOUT



**Antares 846**  
Rigid cervical collar in 2 pieces size  
STANDARD



**Antares 848**  
Rigid cervical collar in 2 pieces size  
MEDIUM



**Antares 842**  
Rigid cervical collar in 2 pieces size  
PEDIATRIC

---