



860 Sirius

"sirius" rigid cervical collar in 1 piece size INFANT



FEATURES

Rigid cervical collar in 1 piece made with transparent polyethylene with E.V.A. washable inner lining fixed by plastic rivets. Adjustable rip-off fastening which different colours indicate the various sizes. Central hole for tracheotomy. Available in 6 different sizes.

MAIN FEATURES

Size	Children
Usage	Immobilization
Trauma	
Pediatric	

SEE ALSO



Antares 854

Cervical collar set 4 pieces with carrying case Art. 842 - 844 - 846 - 848



Sirius 876

Cervical collar set 6 pieces with carrying case