



# 862 Sirius

Rigid cervical collar in 1 piece size PEDIATRIC



## FEATURES

Rigid cervical collar in 1 piece made with transparent polyethylene with E.V.A. washable inner lining fixed by plastic rivets. Adjustable rip-off fastening which different colours indicate the various sizes. Central hole for tracheotomy. Available in 6 different sizes.

## MAIN FEATURES

Size	Children
Usage	Immobilization
Trauma	
Pediatric	

## SEE ALSO



**Antares 854**  
Cervical collar set 4 pieces with carrying case Art. 842 - 844 - 846 - 848



**Sirius 876**  
Cervical collar set 6 pieces with carrying case



**Antares 856**  
Cervical collar set 6 pieces with carrying case



**Sirius 874**  
Cervical collar set 4 pieces with carrying case Art. 862 - 864 - 866 - 868

---

ACCESSORY OF

---



**1366**

Red knapsack "115" with side  
pockets and compartments  
COMPLETE

---