



## 862 Sirius

Rigid cervical collar in 1 piece size PEDIATRIC



### FEATURES

Rigid cervical collar in 1 piece made with transparent polyethylene with E.V.A. washable inner lining fixed by plastic rivets. Adjustable rip-off fastening which different colours indicate the various sizes. Central hole for tracheotomy. Available in 6 different sizes.

### MAIN FEATURES

|           |                |
|-----------|----------------|
| Size      | Children       |
| Usage     | Immobilization |
| Trauma    |                |
| Pediatric |                |

### SEE ALSO



#### Antares 854

Cervical collar set 4 pieces with carrying case Art. 842 - 844 - 846 - 848



#### Sirius 876

Cervical collar set 6 pieces with carrying case



#### Antares 856

Cervical collar set 6 pieces with carrying case



#### Sirius 874

Cervical collar set 4 pieces with carrying case Art. 862 - 864 - 866 - 868

---

ACCESSORY OF

---



**1366**  
Red knapsack "115" with side  
pockets and compartments  
COMPLETE

---