



## 874 Sirius

---

Cervical collar set 4 pieces with carrying case  
Art. 862 - 864 - 866 - 868



---

### FEATURES

Rigid cervical collar in 1 piece made with transparent polyethylene with E.V.A. washable inner lining fixed by plastic rivets. Adjustable rip-off fastening which different colours indicate the various sizes. Central hole for tracheotomy. Available in 6 different sizes.  
Art. 862 - 864 - 866 - 868

---

### MAIN FEATURES

---

Usage

Immobilization

Trauma



SEE ALSO



**Sirius 876**  
Cervical collar set 6 pieces with  
carrying case



**Antares 856**  
Cervical collar set 6 pieces with  
carrying case



**Antares 854**  
Cervical collar set 4 pieces with  
carrying case Art. 842 - 844 - 846 -  
848



**Sirius 862**  
Rigid cervical collar in 1 piece size  
PEDIATRIC



**Sirius 864**  
Rigid cervical collar in 1 piece size  
STOUT



**Sirius 866**  
Rigid cervical collar in 1 piece size  
STANDARD



**Sirius 868**  
Rigid cervical collar in 1 piece size  
MEDIUM