



# 876 Sirius

Cervical collar set 6 pieces with carrying case



## FEATURES

Rigid cervical collar in 1 piece made with transparent polyethylene with E.V.A. washable inner lining fixed by plastic rivets. Adjustable rip-off fastening which different colours indicate the various sizes. Central hole for tracheotomy. Available in 6 different sizes.

## MAIN FEATURES

Usage	Immobilization
Trauma	

SEE ALSO



**Antares 856**  
Cervical collar set 6 pieces with carrying case



**Antares 856**  
Cervical collar set 6 pieces with carrying case



**Sirius 874**  
Cervical collar set 4 pieces with carrying case Art. 862 - 864 - 866 - 868



**Sirius 876**  
Cervical collar set 6 pieces with carrying case



**Sirius 862**  
Rigid cervical collar in 1 piece size PEDIATRIC



**Sirius 864**  
Rigid cervical collar in 1 piece size STOUT



**Sirius 866**  
Rigid cervical collar in 1 piece size STANDARD



**Sirius 868**  
Rigid cervical collar in 1 piece size MEDIUM



**Sirius 860**  
"sirius" rigid cervical collar in 1 piece size INFANT



**Sirius 870**  
Rigid cervical collar in 1 piece size TALL