



880 Halley

Vacuum splint for leg



FEATURES

Produced with absolute reliability materials this splints allow the perfect immobilization of the lower and upper limbs by taking away the inside air. They are supplied with suction valve positioned on the out-side to avoid problems to the patient. It is possible to order the single vacuum splint or the complete set including splints for leg, arm, forearm, little pump for the vacuum and nylon bag. All splints are X-ray transparent.

MAIN FEATURES

Size	Adults
Usage	Immobilization
Trauma	
Length (cm)	109
Width (cm)	67
Material	PVC

IMAGES



SEE ALSO



Halley 886
Vacuum splint set with pump and carrying case



Halley 884
Vacuum splint for forearm



Halley 886
Vacuum splint set with pump and carrying case

ACCESSORIES



889
Little aluminium pump