



882 Halley

Vacuum splint for arm



FEATURES

Produced with absolute reliability materials this splints allow the perfect immobilization of the lower and upper limbs by taking away the inside air. They are supplied with suction valve positioned on the out-side to avoid problems to the patient. It is possible to order the single vacuum splint or the complete set including splints for leg, arm, forearm, little pump for the vacuum and nylon bag. All splints are X-ray transparent.

MAIN FEATURES

Size	Adults
Usage	Immobilization
Trauma	⊌
Length (cm)	75
Width (cm)	61
Material	PVC


IMAGES





SEE ALSO







Halley 880 Vacuum splint for leg



Halley 886Vacuum splint set with pump and carrying case

ACCESSORIES

