



11200 Vital Test


Handheld pulse oximeter complete with adult sensor



FEATURES

Vital Test is simple and professional, it can be used in intensive care units, surgeries, dentist's surgeries, emergency rooms, since it guarantees immediate, accurate and non-invasive monitoring. Vital Test has pleasant and ergonomic shape, it is easy to grip and to be transported. Vital Test allows SpO2 measurement, heart rate, waveform visualization and also last 10 minutes trend value of these two parameters. Supplied with rechargeable battery.

MAIN FEATURES

Cardio	
Width (cm)	7
Height (cm)	15,8
Depth (cm.)	3,7
Weight (kg)	0,25

SEE ALSO



11220 Oxytest
OXYTEST portable pulse oximeter
complete with adult sensor, USB
cable and AAA battery

ACCESSORIES



11204
Support for pulse oxymeter



11208
Pediatric sensor for pulse oximeter
Vital Test



11212
Neonatal sensor for pulse oximeter
Vital Test

STANDARD ACCESSORIES



11206
Adult sensor for pulse oximeter
Vital Test



11210
Handheld pulse oximeter cover