



# 11200 Vital Test


Handheld pulse oximeter complete with adult sensor



## FEATURES

Vital Test is simple and professional, it can be used in intensive care units, surgeries, dentist's surgeries, emergency rooms, since it guarantees immediate, accurate and non-invasive monitoring. Vital Test has pleasant and ergonomic shape, it is easy to grip and to be transported. Vital Test allows SpO2 measurement, heart rate, waveform visualization and also last 10 minutes trend value of these two parameters. Supplied with rechargeable battery.

## MAIN FEATURES

Cardio	
Width (cm)	7
Height (cm)	15,8
Depth (cm.)	3,7
Weight (kg)	0,25

---

SEE ALSO

---



**11220 Oxytest**  
OXYTEST portable pulse oximeter  
complete with adult sensor, USB  
cable and AAA battery

---

ACCESSORIES

---



**11204**  
Support for pulse oximeter



**11208**  
Pediatric sensor for pulse oximeter  
Vital Test



**11212**  
Neonatal sensor for pulse oximeter  
Vital Test

---

STANDARD ACCESSORIES

---



**11206**  
Adult sensor for pulse oximeter  
Vital Test



**11210**  
Handheld pulse oximeter cover