



895 Snake Plus

Independent chambers vacuum mattress with integrated head immobilizer



FEATURES

Mattress with integrated head immobilizer to allow a fast head immobilization. Made of polyurethane treated with fire retardant and in conformity with UNI EN 1865, 14 independent chambers for a perfect and safe immobilization. Usable between - 30°C and + 70°C. Extreme light weight and strenght, 8 ergonomic handles, padded and reinforced, safety belts and a perfect containment of the patient, give this mattress high performances. Usable with X-ray, CT and MRI.

[VIDEO TUTORIAL](#)

MAIN FEATURES

Size	Adults
Usage	Immobilization
Length (cm)	210
Width (cm)	100
Weight (kg)	8,3
Loading capacity (kg)	150
Trauma	
Material	Polyurethane

SEE ALSO



Snake 894
Independent chambers vacuum
mattress

ACCESSORIES



888
Aluminium aspiration pump



890-001
Carrying case for vacuum mattress



896/DF
Abs aspiration pump

STANDARD ACCESSORIES



894-001
Aspiration valve for vacuum splints
and mattress
