



895 Snake Plus

Independent chambers vacuum mattress with integrated head immobilizer




FEATURES

Mattress with integrated head immobilizer to allow a fast head immobilization. Made of polyurethane treated with fire retardant and in conformity with UNI EN 1865, 14 independent chambers for a perfect and safe immobilization. Usable between - 30°C and + 70°C. Extreme light weight and strenght, 8 ergonomic handles, padded and reinforced, safety belts and a perfect containment of the patient, give this mattress high performances. Usable with X-ray, CT and MRI.

[VIDEO TUTORIAL](#)

MAIN FEATURES

| | |
|-----------------------|---|
| Size | Adults |
| Usage | Immobilization |
| Length (cm) | 210 |
| Width (cm) | 100 |
| Weight (kg) | 8,3 |
| Loading capacity (kg) | 150 |
| Trauma |  |
| Material | Polyurethane |

SEE ALSO



Snake 894
Independent chambers vacuum
mattress

ACCESSORIES



888
Aluminium aspiration pump



890-001
Carrying case for vacuum mattress



896/DF
Abs aspiration pump

STANDARD ACCESSORIES



894-001
Aspiration valve for vacuum splints
and mattress