



## 899 Med

Spinal immobilizer



### FEATURES

The spinal immobilizer "MED" is ideal for the immobilization and extraction of a patient blocked in an emergency situation as in case of accident. Thanks to its shape and structure, it allows to fix horizontally and vertically the head and spine to limit all the injuries resulting from a wrong transport of the patient. It is supplied with buckles and belts of different colours for an easy usage, padded cushion and two rip-off straps for forehead and chin. Complete with bag for the transport.

### MAIN FEATURES

Size	Adults
Usage	Immobilization
Trauma	<input checked="" type="checkbox"/>
Width (cm)	85
Height (cm)	83
Bag width (cm)	95
Whole length (cm)	37
Material	Nylon

---

SEE ALSO

---



**Med Ped 9120**  
Pediatric spinal immobilizer

---