



## 9024 Duck

Orange pediatric spinal board



### FEATURES

Designed for immobilization and transport of children with height up to 120 cm. This spinal board has been produced with plastic material at high strenght and studied paying particular attention to the problems connected to daily interventions. It has 4 holes for the quick and total fixing of the head immobilizer and 2 cavities where the board lays on the floor, when base is blocked in traditional way, that allow to avoid damages to the rip-off straps during the usage of the spinal board or its accomodation in the ambulance.

It has 10 handles for the transport.

It allows preliminary disgnostic under X-ray and it's available in 3 different colours: yellow, orange and blue.

It isn't supplied with belts.

### MAIN FEATURES

Size	Children
Usage	Immobilization
Length (cm)	125
Width (cm)	41
Thickness (cm)	5
Weight (kg)	3,7
Loading capacity (kg)	100
Pediatric	<input checked="" type="checkbox"/>
Trauma	<input checked="" type="checkbox"/>
Material	ABS

---

## SEE ALSO

---



**Duck 9022**  
Yellow pediatric spinal board

---

---

## ACCESSORIES

---



**Kids 627**  
Padded head immobilizer –  
pediatric



**Spider Baby 692**  
Immobilization belts for pediatric  
spinal board



**16150-002**  
Storage bag for Ergon stretcher

---