



910 Winner

Self loading stretcher



FEATURES

Article 910 has been designed to offer a practical and safe instrument to the operator in case of emergency. The stretcher, built with anodized aluminium, is strong and handy thanks to its four wheels diam. 150, two fixed and with brakes and two rotating. Important feature of this article is the loading height regulation system. The automatic loading stretcher "WINNER" in fact is supplied with two sturdy perforated plates that allow to change the little front wheels height in order to load the stretcher on ambulances of various model and with different inside fittings.

- Back rest and lying part made with thermoplastic material reinforced to allow cardiac massage. Easy to clean and to disinfect
- Adjustable back rest 0° - 80° moved by a gas piston which allow to fix the back rest safely in any position
- Side protections completely overturnable
- Blocking device for folded legs with a system that avoid the compression of the blocking bar
- Safety lever for the legs positioned near the unlocking device allowing thus the release operation for the loading, keeping the hands on the stretcher
- Vertical legs protected by nylon wedges
- Automatic centring device mounted on rotating wheels. This system automatically block the back wheels in the central position during the loading of the stretcher on the ambulance without having to turn the wheels manually

MAIN FEATURES

Size	Adults
Usage	Transport
Lying surface	Without Trendelenburg
Wheel type	2 Rotating Wheels + 2 Small Fixed Wheels
Type	Self loading stretcher
Height (cm)	87,5
Width (cm)	56
Weight (kg)	34
Loading capacity (kg)	170
Folded height (cm)	26
Max loading height (cm)	65
Length (cm)	194
Min loading height (cm)	45
Material	Anodized aluminium
Length with folded trolley (cm)	162

SEE ALSO



Winner 910 PROOF
Self loading stretcher - certified

ACCESSORIES



Front blocking device
Front blocking device for automatic loading stretcher



Rear fixing system
Rear fixing system for roll-in H=148 mm



950/2
Rear fixing system for roll-in H=160 mm



950/1
Rear fixing system for roll-in H=188 mm



936
Phleboclysis stand for automatic loading stretcher Winner and Frog



939
Instrument support for stretchers



938
Oximeb oxygen cylinder holder



612/AN-2/MEB
Adjustable thoracic abdominal belt orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs



612/N-2/MEB
Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



608/A-MEB
2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)



608/N-MEB
2 pieces black belt with metallic hook and MeBer band (L. 155 cm)



12025
Black multi usage object holder for automatic loading stretcher



7013/A
3 pcs. Skay fire-proof electrowelded mattress



7013/N
3 pcs. Skay fire-proof electrowelded mattress



7014/A
4 pcs. Skay fire-proof electrowelded mattress



12018N
Telescopic foldable phleboclysis stand for self loading stretchers



12070/BG
"COVER OX 2-3 BG"- professional cylinder holder for 2 and 3 lt oxygen bottles



12070/BR
COVER OX2-3 BR professional cylinder holder with belt and orange MeBer logo for 2 and 3 lt oxygen bottles



12072/BG
"COVER OX 5 BG"- professional
cylinder holder for 5 lt oxygen
bottles



12072/BR
"COVER OX 5 BR"- professional
cylinder holder for 5 lt oxygen
bottles



12054
Oximeb Maxi - Oxygen cylinder
holder