



## 910 PROOF Winner

---

Self loading stretcher - certified



---

### FEATURES

Art. 910 PROOF is similar to model 910, the main difference is given by the mean of a rigid steel frame which allows, in coupling with the appropriate fixing system, to absorb those stresses caused by sudden accelerations or decelerations of the rescue vehicle.

Other main features:

- Back rest and lying part made with thermoplastic material reinforced to allow cardiac massage. Easy to clean and to disinfect.
- Adjustable back rest 0° - 80° moved by a gas piston which allows to fix the back rest safely in any position.
- Side protections completely overturnable.
- Blocking device for folded legs with a system that avoids the compression of the blocking bar.
- Safety lever for the legs positioned near the unlocking device allowing thus the release operation for the loading, keeping the hands on the stretcher.
- Vertical legs protected by nylon wedges.
- Automatic centring device mounted on rotating wheels. This system automatically blocks the back wheels in the central position during the loading of the stretcher on the ambulance without having to turn the wheels manually.
- 4 wheels diam. 150, two fixed and with brakes and two rotating.
- Complete with belts and mattress



Device certified according to the European harmonized safety standards UNI EN1865- UNI EN1789.

---

## MAIN FEATURES

---

Size	Adults
Usage	Transport
Lying surface	Without Trendelenburg
Wheel type	2 Rotating Wheels + 2 Small Fixed Wheels
Type	
10G Certification	
Height (cm)	87,5
Width (cm)	56
Weight (kg)	37
Loading capacity (kg)	170
Folded height (cm)	26
Max loading height (cm)	65
Min loading height (cm)	62
Length (cm)	194
Length with folded trolley (cm)	155
Material	Anodized aluminium
Filter feature	

## ACCESSORIES



**955 PROOF**  
EN 1789 fixing system for Proof roll-in stretchers



**936**  
Phleboclysis stand for automatic loading stretcher Winner and Frog



**939**  
Instrument support for stretchers



**938**  
Oximeb oxygen cylinder holder



**12025**  
Black multi usage object holder for automatic loading stretcher



**7013/A**  
3 pcs. Skay fire-proof electrowelded mattress



**7013/N**  
3 pcs. Skay fire-proof electrowelded mattress



**7014/A**  
4 pcs. Skay fire-proof electrowelded mattress



**12018N**  
Telescopic foldable phleboclysis stand for self loading stretchers



**612/AN-2/MEB**  
Adjustable thoracic abdominal belt orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs



**612/N-2/MEB**  
Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



**608/A-MEB**  
2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)



**608/N-MEB**  
2 pieces black belt with metallic hook and MeBer band (L. 155 cm)



**12070/BG**  
"COVER OX 2-3 BG"- professional cylinder holder for 2 and 3 lt oxygen bottles



**12070/BR**  
COVER OX2-3 BR professional cylinder holder with belt and orange MeBer logo for 2 and 3 lt oxygen bottles



**12072/BG**  
"COVER OX 5 BG"- professional cylinder holder for 5 lt oxygen bottles



**12072/BR**

"COVER OX 5 BR"- professional cylinder holder for 5 lt oxygen bottles

---

## STANDARD ACCESSORIES

---



**7013/A**

3 pcs. Skay fire-proof electrowelded mattress



**608/A-MEB**

2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)



**608/N-MEB**

2 pieces black belt with metallic hook and MeBer band (L. 155 cm)



**612/AN-2/MEB**

Adjustable thoracic abdominal belt orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs



**612/N-2/MEB**

Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs

---