



## 9100 Bunny

Pediatric vacuum mattress with independent chambers



### FEATURES

Mattress made of polyurethane with fire retardant treatment, 14 independent chambers for a perfect and safe immobilization. Usable between -30°C and +70°C. Extreme light weight and strong, 8 ergonomic handles padded and reinforced, safety belts and a perfect containment of the patient, give this mattress high performances. Usable with X-ray, CR and MRI. Weight: 3 kg.

### MAIN FEATURES

Size	Children
Usage	Immobilization
Length (cm)	120
Width (cm)	70
Weight (kg)	3
Pediatric	<input checked="" type="checkbox"/>
Trauma	<input checked="" type="checkbox"/>
Material	Polyurethane

### SEE ALSO



#### Bunny Plus 9102

Independent chambers vacuum mattress pediatric with integrated head immobilizer

---

## ACCESSORIES

---



**888**  
Aluminium aspiration pump



**896/DF**  
Abs aspiration pump



**9100-01**  
Carrying bag for pediatric vacuum mattress

---

## STANDARD ACCESSORIES

---



**894-001**  
Aspiration valve for vacuum splints and mattress

---