



## 9102 Bunny Plus

Independent chambers vacuum mattress pediatric with integrated head immobilizer



### FEATURES

Mattress with integrated head immobilizer to allow a fast head immobilization. Made of polyurethane treated with fire retardant, **14 independent chambers** for a perfect and safe immobilization. Usable between - 30°C and + 70°C. Extreme light weight and strenght, 8 ergonomic handles padded and reinforced, safety belts and a perfect containment of the patient, give this mattress high performances. Usable with X-ray, CT and MRI.

### MAIN FEATURES

Size	Children
Usage	Immobilization
Length (cm)	120
Width (cm)	70
Weight (kg)	3,5
Pediatric	<input checked="" type="checkbox"/>
Trauma	<input checked="" type="checkbox"/>
Material	Polyurethane

---

## SEE ALSO

---



**Bunny 9100**  
Pediatric vacuum mattress with  
independent chambers

---

---

## ACCESSORIES

---



**888**  
Aluminium aspiration pump



**896/DF**  
Abs aspiration pump

---

---

## STANDARD ACCESSORIES

---



**894-001**  
Aspiration valve for vacuum splints  
and mattress

---