



## 914 PROOF Winner

---

Certified self loading stretcher with Trendelenburg and Fowler



---

### FEATURES

Art. 914 PROOF is similar to model 914, the main difference is given by the mean of a rigid steel frame which allows, in coupling with the appropriate fixing system, to absorb those stresses caused by sudden accelerations or decelerations of the rescue vehicle. Complete with mattress and belts.

This article has been designed to offer a practical and safe instrument to the operator in case of emergency. The stretcher, built with anodized aluminium, is strong and handy thanks to its four wheels diam. 150, two fixed with brakes and two rotating. Important feature of this article is the loading height regulation system. This stretcher in fact is supplied with two sturdy perforated plates that allow to change the little front wheels height in order to load the stretcher on ambulances of various model and with different inside fittings.

It also has the possibility of having Trendelenburg and Fowler positions becoming this way a more complete and functional support for every emergency situation. Back rest and lying part made with thermoplastic material reinforced to allow cardiac massage.

Easy to clean and to disinfect. Adjustable back rest 0° - 80°, moved by a gas piston that allow to fix the back rest safely in any position. Trendelenburg and Fowler positions moved by gas pistons. Side protections completely overturnable. Blocking device for folded legs with a system that avoid the compression of the blocking bar. Safety lever for the legs positioned near the unlocking device allowing thus the release operation for the loading, keeping the hands on the stretcher. Vertical legs protected by nylon wedges. Self centring device mounted on rotating wheels.

This system automatically blocks the back wheels in the central position during the loading of the stretcher on the ambulance avoiding to turn the wheels manually.




Device certified according to the European harmonized safety standards UNI EN1865 - UNI EN 1789.

---

## MAIN FEATURES

---

|                                 |   |
|---------------------------------|---|
| Usage                           | Transport   |
| Lying surface                   | With Trendelenburg and Fowler   |
| Wheel type                      | 2 Rotating Wheels + 2 Small Fixed Wheels  |
| Type                            | Self loading stretcher  |
| 10G Certification               |  |
| Height (cm)                     | 87,5  |
| Width (cm)                      | 56  |
| Weight (kg)                     | 40  |
| Loading capacity (kg)           | 170   |
| Folded height (cm)              | 26  |
| Max loading height (cm)         | 65  |
| Min loading height (cm)         | 62  |
| Length (cm)                     | 194   |
| Length with folded trolley (cm) | 155   |
| Material                        | Anodized aluminium  |

---

## SEE ALSO

---



**Winner 910 PROOF**  
Self loading stretcher - certified

---

ACCESSORIES



**7014/A**  
4 pcs. Skay fire-proof electrowelded  
mattress



**955 PROOF**  
EN 1789 fixing system for Proof roll-  
in stretchers



**936**  
Phlebotomy stand for automatic  
loading stretcher Winner and Frog



**939**  
Instrument support for stretchers



**938**  
Oximeb oxygen cylinder holder



**612/N-2/MEB**  
Adjustable thoracic abdominal belt  
black with MeBer band with 4  
points fastening for self loading  
stretchers and sedan chairs



**608/A-MEB**  
2 pieces orange belt with metallic  
hook and MeBer band (L. 155 cm)



**608/N-MEB**  
2 pieces black belt with metallic  
hook and MeBer band (L. 155 cm)



**12025**  
Black multi usage object holder for  
automatic loading stretcher



**612/AN-2/MEB**  
Adjustable thoracic abdominal belt  
orange/black with MeBer band and  
4 points fastening for self loading  
stretchers and sedan chairs



**7013/A**  
3 pcs. Skay fire-proof electrowelded  
mattress



**7013/N**  
3 pcs. Skay fire-proof electrowelded  
mattress



**12018N**  
Telescopic foldable phlebotomy  
stand for self loading stretchers



**12070/BG**  
"COVER OX 2-3 BG"- professional  
cylinder holder for 2 and 3 lt oxygen  
bottles



**12070/BR**  
COVER OX2-3 BR professional  
cylinder holder with belt and  
orange MeBer logo for 2 and 3 lt  
oxygen bottles



**12072/BG**  
"COVER OX 5 BG"- professional  
cylinder holder for 5 lt oxygen  
bottles



**12072/BR**  
"COVER OX 5 BR"- professional  
cylinder holder for 5 lt oxygen  
bottles



**12054**  
Oximeb Maxi - Oxygen cylinder  
holder

---

## STANDARD ACCESSORIES

---



**608/A-MEB**  
2 pieces orange belt with metallic  
hook and MeBer band (L. 155 cm)



**608/N-MEB**  
2 pieces black belt with metallic  
hook and MeBer band (L. 155 cm)



**612/AN-2/MEB**  
Adjustable thoracic abdominal belt  
orange/black with MeBer band and  
4 points fastening for self loading  
stretchers and sedan chairs



**612/N-2/MEB**  
Adjustable thoracic abdominal belt  
black with MeBer band with 4  
points fastening for self loading  
stretchers and sedan chairs



**7013/A**  
3 pcs. Skay fire-proof electrowelded  
mattress

---