



9150 X-Lock

Ergonomic head immobilizer for Ergon stretcher



FEATURES

One piece head immobilizer with integrated head strap, for a fast head immobilization directly on the stretcher. Made with laminated polyesters and bi-density polyurethane foams.

MAIN FEATURES

Size	Adults
Usage	Immobilization
Trauma	
Length (cm)	43
Width (cm)	27
Height (cm)	140

SEE ALSO



Ergon 16150/V
Variable geometry comfort
stretcher - green/blue



Ergon 16150/G
Variable geometry comfort
stretcher - yellow/blue

ACCESSORIES



Art. 9150-001
Blue carrying case for X-Lock
ergonomic head immobilizer
