



936/A

Phleboclysis stand for stretchers Art. 1510 Proof and Art. 1515 Proof



FEATURES

Telescopic phleboclysis stand for stretcher. Thanks to its fixing system it is possible to fix it in every point on the stretcher's frame. **Dimensions: minimum height: cm 47,5 maximum height: cm 130,5**

MAIN FEATURES

Min height (cm)	47,5
Max height (cm)	130,5
Weight (kg)	0,5
Material	Aluminium alloy

IMAGES



SEE ALSO



936

Phlebotomy stand for automatic loading stretcher Winner and Frog



12018N

Telescopic foldable phlebotomy stand for self loading stretchers
