



938

Oximeb oxygen cylinder holder



FEATURES

Cylinder holder produced with strong, antitear nylon and inside padding to prevent damages from knocks. It is supplied with a handle for the transport and a side pocket for tube and mask.

Can be used with cylinders of 2 liters and, thanks to a handy fastening, it can be mounted on the side protection of the automatic loading stretchers and on the back rest of the sedan chairs.

MAIN FEATURES

Weight (kg)	0,3
Length (cm)	31
Width (cm)	19
Depth (cm.)	13
Material	Nylon

SEE ALSO



12054
Oximeb Maxi - Oxygen cylinder
holder

ACCESSORY OF



Mercury Lite 7082/4RG Proof
Certified aluminium variable height
stretcher high version with 4
rotating wheels



Mercury Cinque 7075/4RG PROOF
Certified 5 levels self-loading
stretcher with 4 swivel wheels



**Mercury Lite Cinque 7095/4RG
Proof**
Certified aluminium 5 levels self-
loading stretcher with 4 swivel
wheels



**Mercury Lite Cinque 7097/4RG
Proof**
Certified aluminium 5 levels self-
loading stretcher with 4 swivel
wheels - high version



Mercury EL 7092 PROOF
Certified self-loading stretcher



7265 Proof
Aluminium self loading stretcher
"Frog lite Advance" with
extractable handles and
Trendelenburg