



939

Instrument support for stretchers



FEATURES

Instrument support to be fixed on Me.Ber. stretcher frames¹ to transport medical equipments, defibrillators and monitors. It is supplied with belts to fix devices and is foldable to be stored inside ambulances.

Dimensions: open cm 34x56x37h, closed height cm 8

MAIN FEATURES

Width (cm)	56
Depth (cm.)	34
Height (cm)	37
Epaisseur pliée (cm)	8
Weight (kg)	5,6
Loading capacity (kg)	30

IMAGES



SEE ALSO



Mercury Lite 7082/4RG Proof
Certified aluminium variable height
stretcher high version with 4
rotating wheels

ACCESSORY OF



Mercury Cinque 7075/4RG PROOF
Certified 5 levels self-loading
stretcher with 4 swivel wheels



Mercury Cinque 7077/4RG PROOF
Certified 5 levels self-loading
stretcher with 4 swivel wheels -
high version



**Mercury Lite Cinque 7095/4RG
Proof**
Certified aluminium 5 levels self-
loading stretcher with 4 swivel
wheels



**Mercury Lite Cinque 7097/4RG
Proof**
Certified aluminium 5 levels self-
loading stretcher with 4 swivel
wheels - high version



Mercury EL 7092 PROOF
Certified self-loading stretcher



Mercury EL 7094 PROOF
Certified self-loading stretcher -
high version



7265 Proof
Aluminium self loading stretcher
"Frog lite Advance" with
extractable handles and
Trendelenburg