



950/4 Rear fixing system

Rear fastening system for ambulance cots H=129 mm



FEATURES

Built in steel, this system are installed on the ambulance floors to safely lock an ambulance cot. It is available in 4 different heights.

MAIN FEATURES

Height (cm)

12,9

SEE ALSO



950/1

Rear fixing system for roll-in H=188 mm



950/2

Rear fixing system for roll-in H=160 mm



Rear fixing system

Rear fixing system for roll-in H=148 mm

ACCESSORY OF



Rescue stretcher 638

Collapsible rescue stretcher with adjustable backrest



Rescue stretcher 640

Collapsible rescue stretcher
