



## 955 PROOF

---

EN 1789 fixing system for Proof roll-in stretchers



---

### FEATURES

Designed to guarantee maximum stability of the stretcher on board the rescue vehicle, they are made entirely of stainless steel for longer life and corrosion resistance. Their easy-to-install structure ensures that the stretcher or trolley, and patient are securely restrained in all directions. Art. 955 Proof is suitable for fixing all our 10G roll-in stretchers whilst art. 956 Proof is suitable for fixing "Titanus" trolleys.



Device certified according to the European harmonized standards UNI EN 1789.

Height	Width	Length	Weight
13 cm	48 cm	177 cm	7,7 kg

---

### MAIN FEATURES

Height (cm)	13
Length (cm)	177
Weight (kg)	7,7

---

## SEE ALSO

---



**956 PROOF**  
Fixing system for Titanus automatic  
loading trolley

---

---

## ACCESSORIES

---



**955-001**  
Front reverse-plate for 955-956  
Proof certified fixing system



**955-002**  
Rear reverse-plate for 955-956 Proof  
certified fixing system

---

---

## ACCESSORY OF

---



**MBS-01 MEBER BARIATRIC STRETCHER 7270**  
Self loading stretcher for bariatric transport



**Mercury Cinque 7075/4RG PROOF**  
Certified 5 levels self-loading stretcher with 4 swivel wheels



**Mercury Cinque 7077/4RG PROOF**  
Certified 5 levels self-loading stretcher with 4 swivel wheels - high version



**Mercury Lite 7080/4RG Proof**  
Ambulance cot in aluminium with variable height and four swivel castors



**Mercury Lite 7082/4RG Proof**  
Certified aluminium variable height stretcher high version with 4 rotating wheels



**Mercury 7070/4RG PROOF**  
Stainless steel variable height stretcher certified with 4 rotating wheels



**Mercury 7072/4RG PROOF**  
Certified stainless steel variable height stretcher high version with 4 rotating wheels



**Frog 7240 PROOF**  
"FROG PLUS" self loading stretcher with trendelenburg certified



**Frog Lite 7260 PROOF**  
"FROG LITE" self loading aluminium stretcher with trendelenburg certified



**Frog 7210 PROOF**  
"FROG" yellow self loading stretcher with trendelenburg certified



**Mercury Lite Cinque 7095/4RG Proof**  
Certified aluminium 5 levels self-loading stretcher with 4 swivel wheels



**Mercury Lite Cinque 7097/4RG Proof**  
Certified aluminium 5 levels self-loading stretcher with 4 swivel wheels - high version



**Mercury Lite 7082/4RG Proof**  
Certified aluminium variable height stretcher high version with 4 rotating wheels



**Mercury EL 7094 PROOF**  
Certified self-loading stretcher - high version



**7265 Proof**  
Aluminium self loading stretcher "Frog lite Advance" with extractable handles and Trendelenburg