



956 PROOF

Fixing system for Titanus automatic loading trolley



FEATURES

Designed to guarantee maximum stability of the stretcher on board the rescue vehicle, they are made entirely of stainless steel for longer life and corrosion resistance. Their easy-to-install structure ensures that the stretcher or trolley, and patient are securely restrained in all directions. Art. 955 Proof is suitable for fixing all our 10G roll-in stretchers whilst art. 956 Proof is suitable for fixing "Titanus" trolleys.



Device certified according to the European harmonized standards UNI EN 1789.

Height	Width	Length	Weight
13 cm	42 cm	175 cm	8,5 kg

MAIN FEATURES

Length (cm)	175
Width (cm)	42
Height (cm)	13
Weight (kg)	8,5

SEE ALSO



955 PROOF
EN 1789 fixing system for Proof roll-in stretchers

ACCESSORIES



955-001
Front reverse-plate for 955-956 Proof certified fixing system



955-002
Rear reverse-plate for 955-956 Proof certified fixing system

ACCESSORY OF



Titanus 7060/4RG PROOF
Certified stainless steel variable height trolley with four rotating wheels



Titanus 7062/4RG PROOF
Stainless steel variable height trolley with four rotating wheels high version