



## 965 Front blocking device

---

Front blocking device for automatic loading stretcher



---

### FEATURES

Blocking device to be fixed on the ambulance platform that allows to guide and block the automatic loading stretcher in the front part. Made of stainless steel.

---

### MAIN FEATURES

Length (cm)	51,5
Width (cm)	48
Height (cm)	10
Weight (kg)	1,4

---

### SEE ALSO



**Winner 914**  
Self loading stretcher with  
trendelenburg and fowler



**Winner 910**  
Self loading stretcher