



959

Wall mount system for variable geometry stretcher Ergon



FEATURES

Wall mount system for Ergon scoop stretcher, made in thermoplastics. This device ensures a safe and quick installation on any emergency vehicle, vertically.

Device certified according to the European harmonized standards UNI EN 1789



MAIN FEATURES

Length (cm)	51
Height (cm)	2,5
Depth (cm.)	9
Weight (kg)	1

IMAGES



SEE ALSO



Ergon 16150/G
Variable geometry comfort
stretcher - yellow/blue



Ergon 16150/V
Variable geometry comfort
stretcher - green/blue
