



959

Wall mount system for variable geometry stretcher Ergon



FEATURES

Wall mount system for Ergon scoop stretcher, made in thermoplastics. This device ensures a safe and quick installation on any emergency vehicle, vertically.

Device certified according to the European harmonized standards UNI EN 1789



MAIN FEATURES

| | |
|-------------|-----|
| Length (cm) | 51 |
| Height (cm) | 2,5 |
| Depth (cm.) | 9 |
| Weight (kg) | 1 |


IMAGES



SEE ALSO



Ergon 16150/G
Variable geometry comfort
stretcher - yellow/blue



Ergon 16150/V
Variable geometry comfort
stretcher - green/blue
