



## 960

---

“Extra fixing system” - fixing system for chair Extra



---

### FEATURES

Fixing device designed to ensure a permanent fixation of the Extra chair inside the ambulance patient compartment. It is made of corrosion resistant and easy to clean materials (stainless and thermoplastic). Extra fixing system enables to lock the whole Extra range in a simple, fast and effective way. The operator can easily lock and unlock the chair, thanks to the innovative and exclusive monobloc configuration.

Device certified according to the European harmonized standards UNI EN 1789.



---

### MAIN FEATURES

Width (cm)	23
Height (cm)	42
Depth (cm.)	20
Weight (kg)	5,8

---

## IMAGES

---



---

## SEE ALSO

---



**975**

"Ergolift fixing system" - Fixing system for chair Extra Ergolift

---

## ACCESSORY OF

---



**Extra 670/BR**

Track assisted stair chair with arms



**679**

Narrow foldable sedan chair



**Extra Lite Standard 676**

Foldable carrying chair with 4 wheels, footrest and arms



**678 Extra Lite Basic**

Extra Lite Basic, foldable carrying chair

---