



16160/G Ergon Junior

Pediatric variable geometry comfort stretcher - yellow/blue



FEATURES

[Visit Ergon website](#)

The greatest product for the smallest patients

ERGON shares its advanced technology with the junior model. All benefits and main features of the adult stretcher are preserved in its Junior version.

The **Flexilock locking system** was moved from the centre to the side, in order to allow operations without invading the head space. The large blade allows for maximum comfort even to pediatric patients.



MAIN FEATURES


| | |
|-----------------------|--|
| Size | Children |
| Usage | Rescue Immobilization Transport |
| Height (cm) | 7,5 |
| Width (cm) | 43 |
| Length (cm) | 129 |
| Max length (cm) | 153 |
| Folded length (cm) | 85 |
| Epaisseur pliée (cm) | 9 |
| Weight (kg) | 8,4 |
| Loading capacity (kg) | 170 |
| Type | |
| Pediatric | <input checked="" type="checkbox"/> |
| Trauma | <input checked="" type="checkbox"/> |
| Alpine | <input checked="" type="checkbox"/> |
| Material | HC3 - High Comfort Composite Cell Aluminium alloy |

SEE ALSO



Ergon Junior 16160/V
Pediatric variable geometry
comfort stretcher - green/blue

ACCESSORIES



X-FIX 9156
Restraint system with straps and
carabiners



959
Wall mount system for variable
geometry stretcher Ergon



16150-002
Storage bag for Ergon stretcher



16150-001
Winch bag for Ergon

STANDARD ACCESSORIES



584
Black bag with 3 orange straps in 2
pcs and quick-release plastic
buckle for Ergon stretcher
