



16150-001

Winch bag for Ergon



FEATURES

Designed to facilitate rescue operations, this case for the Ergon scoop stretcher ensures maximum operativity in hostile environment. Equipped with 4 shoulder-straps. The Ergon stretcher (disassembled) fits with all of its accessories inside of the case, strongly reducing its volume. The EN566 hook in the centre ensures safe winch operation. Size: **cm 124x32x16**

MAIN FEATURES

Length (cm)	124
Width (cm)	32
Loading capacity (kg)	16

SEE ALSO



16150-002
Storage bag for Ergon stretcher

ACCESSORY OF



Ergon 16150/G
Variable geometry comfort
stretcher - yellow/blue



Ergon 16150/V
Variable geometry comfort
stretcher - green/blue
