



## 16150-001

---

Winch bag for Ergon



---

### FEATURES

Designed to facilitate rescue operations, this case for the Ergon scoop stretcher ensures maximum operativity in hostile environment. Equipped with 4 shoulder-straps. The Ergon stretcher (disassembled) fits with all of its accessories inside of the case, strongly reducing its volume. The EN566 hook in the centre ensures safe winch operation. Size: **cm 124x32x16**

---

### MAIN FEATURES

Length (cm)	124
Width (cm)	32
Loading capacity (kg)	16

---

### SEE ALSO



**16150-002**  
Storage bag for Ergon stretcher

---

---

## ACCESSORY OF

---



**Ergon 16150/G**  
Variable geometry comfort  
stretcher - yellow/blue



**Ergon 16150/V**  
Variable geometry comfort  
stretcher - green/blue

---