



# 16150-002

Storage bag for Ergon stretcher



## FEATURES

Measuring as little as 11 cm in depth, this bag is ideal for Ergon storage, together with its accessories. It was designed for storage purposes or equipped on emergency vehicles. Once opened on the ground, the bag works as an insulating layer, allowing to safely prepare medicines or operate electro-medical devices. Size: **cm 122x47x8**

## MAIN FEATURES

Length (cm)	122
Width (cm)	47
Depth (cm.)	8

## IMAGES



---

SEE ALSO

---



**16150-001**  
Winch bag for Ergon

---

ACCESSORY OF

---



**Ergon 16150/G**  
Variable geometry comfort  
stretcher - yellow/blue



**Ergon 16150/V**  
Variable geometry comfort  
stretcher - green/blue

---