



7122/AZ Gecko

“GECKO” self loading stretcher entry-level light blue



FEATURES

Studied to optimize quality-price ratio. Build with painted stainless steel. This device grants all features of other self loading stretchers as: adjustable back rest, 2 fixed wheels with brake and 2 rotating wheels Ø 200 mm self centring, folding side protections. Supplied with belts.

MAIN FEATURES

Type	
Usage	Transport
Lying surface	Without Trendelenburg
Height (cm)	90,5
Max loading height (cm)	62
Min loading height (cm)	62
Folded height (cm)	30
Width (cm)	56
Length (cm)	197
Weight (kg)	40,2
Loading capacity (kg)	170
Material	Painted steel

IMAGES



SEE ALSO



Gecko 7122/V
"GECKO" self loading stretcher
entry-level green

ACCESSORIES



939
Instrument support for stretchers



Rear fixing system
Rear fixing system for roll-in H=148 mm



Front blocking device
Front blocking device for automatic loading stretcher



938
Oximeb oxygen cylinder holder



7016/A
PVC orange electrowelded mattress (one-piece)



612/AN-2/MEB
Adjustable thoracic abdominal belt orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs



612/N-2/MEB
Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



12080-001
Monting kit for bariatric platform on self loading stretcher



12010
Phleboclysis stand for stretchers Mercury and Gecko



12025
Black multi usage object holder for automatic loading stretcher



608/A-MEB
2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)



608/N-MEB
2 pieces black belt with metallic hook and MeBer band (L. 155 cm)



12080
"MPB" wide bariatric platform for self-loading stretchers Frog Mercury and Gecko complete of belts with reel and fixing for stretcher

STANDARD ACCESSORIES



603/MEB
Nylon bag containing 3 belts Art. 607/Ameb