



7122/V Gecko

"GECKO" self loading stretcher entry-level green



FEATURES

Studied to optimize quality-price ratio. Build with painted stainless steel. This device grants all features of other self loading stretchers as: adjustable back rest, 2 fixed wheels with brake and 2 rotating wheels Ø 200 mm self centring, folding side protections. Suplied with belts.

MAIN FEATURES

Туре	
Usage	Transport
Lying surface	Without Trendelenburg
Height (cm)	90,5
Max loading height (cm)	62
Min loading height (cm)	62
Folded height (cm)	30
Width (cm)	56
Length (cm)	197
Weight (kg)	37
Loading capacity (kg)	170
Material	Painted steel

IMAGES





SEE ALSO



Gecko 7122/AZ"GECKO" self loading stretcher entry-level light blue

ACCESSORIES



7016/APVC orange electrowelded mattress (one-piece)



939 Instrument support for stretchers



Rear fixing system
Rear fixing system for roll-in H=148
mm



Front blocking device Front blocking device for automatic loading stretcher



938 Oximeb oxygen cylinder holder



612/AN-2/MEB

Adjustable thoracic abdominal belt orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs



612/N-2/MEBAdjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



12080-001 Monting kit for bariatric platform on self loading stretcher



12025
Black multi usage object holder for automatic loading stretcher



12080
"MPB" wide bariatric platform for self-loading stretchers Frog Mercury and Gecko complete of belts with reel and fixing for stretcher



12010Phleboclysis stand for stretchers
Mercury and Gecko



608/A-MEB2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)



2 pieces black belt with metalic hook and MeBer band (L. 155 cm)



12070/BG"COVER OX 2-3 BG"- professional cylinder holder for 2 and 3 lt oxygen bottles



12070/BR COVER OX2-3 BR professional cylinder holder with belt and orange MeBer logo for 2 and 3 lt oxygen bottles



12072/BG"COVER OX 5 BG"- professional cylinder holder for 5 lt oxygen bottles



12072/BR
"COVER OX 5 BR"- professional
cylinder holder for 5 lt oxygen
bottles

STANDARD ACCESSORIES



603/MEBNylon bag containing 3 belts Art.
607/Ameb