



7122/V Gecko

“GECKO” self loading stretcher entry-level green



FEATURES

Studied to optimize quality-price ratio. Build with painted stainless steel. This device grants all features of other self loading stretchers as: adjustable back rest, 2 fixed wheels with brake and 2 rotating wheels Ø 200 mm self centring, folding side protections. Suplied with belts.

MAIN FEATURES

Type	
Usage	Transport
Lying surface	Without Trendelenburg
Height (cm)	90,5
Max loading height (cm)	62
Min loading height (cm)	62
Folded height (cm)	30
Width (cm)	56
Length (cm)	197
Weight (kg)	40,2
Loading capacity (kg)	170
Material	Painted steel

IMAGES



SEE ALSO



Gecko 7122/AZ
"GECKO" self loading stretcher
entry-level light blue

ACCESSORIES



7016/A
PVC orange electrowelded
mattress (one-piece)



939
Instrument support for stretchers



Rear fixing system
Rear fixing system for roll-in H=148
mm



Front blocking device
Front blocking device for
automatic loading stretcher



938
Oximeb oxygen cylinder holder



612/AN-2/MEB
Adjustable thoracic abdominal belt
orange/black with MeBer band and
4 points fastening for self loading
stretchers and sedan chairs



612/N-2/MEB
Adjustable thoracic abdominal belt
black with MeBer band with 4
points fastening for self loading
stretchers and sedan chairs



12080-001
Monting kit for bariatric platform
on self loading stretcher



12025
Black multi usage object holder for
automatic loading stretcher



12080
"MPB" wide bariatric platform for
self-loading stretchers Frog
Mercury and Gecko complete of
belts with reel and fixing for
stretcher



12010
Phleboclysis stand for stretchers
Mercury and Gecko



608/A-MEB
2 pieces orange belt with metallic
hook and MeBer band (L. 155 cm)



608/N-MEB
2 pieces black belt with metallic
hook and MeBer band (L. 155 cm)



12070/BG
"COVER OX 2-3 BG"- professional
cylinder holder for 2 and 3 lt oxygen
bottles



12070/BR
COVER OX2-3 BR professional
cylinder holder with belt and
orange MeBer logo for 2 and 3 lt
oxygen bottles



12072/BG
"COVER OX 5 BG"- professional
cylinder holder for 5 lt oxygen
bottles



12072/BR

"COVER OX 5 BR"- professional
cylinder holder for 5 lt oxygen
bottles

STANDARD ACCESSORIES



603/MEB

Nylon bag containing 3 belts Art.
607/Ameb
