



## 891 "Flake"


Independent chambers vacuum mattress



### FEATURES

Vacuum mattress with 3 independent chambers for an excellent and fast immobilization. Extremely light and sturdy, 8 handles, restraining belts and a perfect patient encirclement, give to the mattress unique solutions. Unmatchable vacuum durability without any deflection. X-ray transparent mattress.

### MAIN FEATURES

Size	Adults
Usage	Immobilization
Length (cm)	200
Width (cm)	90
Weight (kg)	6,7
Loading capacity (kg)	150
Trauma	
Material	PVC-coated polyester

### SEE ALSO



**891/SG**  
Flake basic independent chambers  
vacuum mattress



**Chesterfield 9114**  
Vacuum mattress with  
independent chambers with  
central welding

---

## ACCESSORIES

---



**888**  
Aluminium aspiration pump



**890-001**  
Carrying case for vacuum mattress



**896/DF**  
Abs aspiration pump

---

## STANDARD ACCESSORIES

---



**894-001**  
Aspiration valve for vacuum splints  
and mattress